

## Instructions:

- If freezing the milk, remove excess air from the bag and seal along the top.
- Be careful to only fill 6 fl oz - overfilled bags may leak during freezing.
- Lay the filled bag flat in the freezer with nothing on top.
- 

## Warning:

- Milk bags are single use only.
- Do not microwave a milk bag.
- It is not recommended to microwave formula or breast milk.
- Keep milk bags and packaging material away from children and out of their reach.

## Instrucciones:

- Si congela la leche, elimine el exceso de aire de la bolsa y selle a lo largo de la parte superior.
- Tenga cuidado de solo llenar 6 onzas líquidas, porque las bolsas muy llenas podrían derramarse durante la congelación.
- Coloque la bolsa con la leche en el congelador sin nada encima.

## Advertencia:

- Las bolsas de leche son desechables.
- No caliente la bolsa de leche en el microondas.
- No se recomienda calentar la fórmula o la leche materna en el microondas.
- Mantenga las bolsas de leche y los empaques fuera del alcance de los niños.



Store breast milk or formula in the fridge or freezer.

Worry free, leak-resistant, double seal.

Each 6 fl oz bag is BPA free and pre-sterilized.

Available in 30 and 100 packs